



Benefits

- ✓ Provide unambiguous judgment of the test results
- ✓ Digital record keeping for evidential use

Features

- Android 7.1
- Touch Screen
- 100,000 tests records storage
- Thermal Printer
- 2D Barcode Scanner
- Wifi, 4G, Bluetooth, NFC
- GPS
- Maximum 9 drugs detected in one cassette
- Saliva, Hair, Urine & Surface Residue Test

Specification

CPU	Dual-Core 1.8G & Quad-Core 1.4GHz
Display	5.5" Touch Screen (720x1280)
Operating System	Android 7.1
Dimension	Handset: 21.65 x 8.65 x 5.65cm, With Docking: 25 x 22 x 12.2cm
Weight	Handset: 562g (with thermal printer), With Docking: 666g
Battery	7.4V / 3,000mAh Lithium Battery
Battery Life	580 tests (with printing), 800 tests (without printing)
Camera	13 Mega Pixel
RAM	DDR3 2GB
Storage	EMMC 16GB, 100,000 tests records
Working Temperature	0°C to 50°C
Storage Temperature	-20°C to 70°C
Barcode	3095 2D Barcode scanner
GPS	Support
Printer	Built-In Thermal Printer
Internet	4G, NFC(Optional), WIFI
Max. no. of drugs detected	9 (Saliva/Hair), 14 (Urine), 6 (Surface)



C4 Development Ltd

Rm. 1306, 13/F Winful Centre, 30 Shing Yip Street, K.T., Hong Kong

☎ 852-34210580 📠 852-34210583 ✉ cshk@c4.com.hk

🌐 <https://www.aquilascan.com/>

Drug Test Analyzer

WDTP-80



Drug Test Cassette

Urine Test



Storage Temperature:
4°C - 30°C
Expiry Date:
24 months


Hair Test



Storage Temperature:
2°C - 30°C
Expiry Date:
24 months


Drug Test Cassette

Saliva Test



Storage
Temperature:
2°C - 30°C
Expiry Date:
18 to 24 months

Surface Residue Test



Storage
Temperature:
4°C - 30°C
Expiry Date:
18 months

Test Results

Test Results

T1A THC No drug detected
 T2A OPI No drug detected
 T3A AMP No drug detected
 T1B COC No drug detected
 T2B THC No drug detected
 T3B MET No drug detected
 T2C KET No drug detected

Test Results

T1A THC No drug detected
 T2A OPI Lab test required
 T3A AMP No drug detected
 T1B COC No drug detected
 T2B THC No drug detected
 T3B MET No drug detected
 T2C KET No drug detected

Available Drugs

Drug	
AMP	Amphetamines
THC	Marijuana
COC	Cocaine
PCP	Phencyclidine
MOR	Morphine
MET	Methamphetamine
MTD	Methadone
BAR	Barbiturate
BZO	Benzodiazepines
TCA	Nortriptyline
MDMA	Ecstasy
BUO	Buprenorphine

Drug	
BUP	Buprenorphine
EDDP	Methadone Metabolite
PPX	Prophylaxis
ETG	Ethylglucoronide
K2	Synthetic Cannabis
TML	Tramadol
MQL	Methaqualone
COT	Cotinine
FEN	Fentanyl
OXY	Oxycodone
KET	Ketamine
OPI	Opiates

Note: Customized drug panels available